

Recombinant Human HEBP2

Catalog No: #GP13285

Package Size: #GP13285-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human HEBP2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 6-205 amino acids of human HEBP2
Target Name	heme binding protein 2
Other Names	PP23; SOUL; C6orf34; C6ORF34B
Accession No.	Swissprot:Q9Y5Z4Gene Accession:BC008205
Uniprot	Q9Y5Z4
GeneID	23593;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene is found predominately in the cytoplasm, where it plays a role in the collapse of mitochondrial membrane potential (MMP) prior to necrotic cell death. The encoded protein enhances outer and inner mitochondrial membrane permeabilization, especially under conditions of oxidative stress.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only