

SDCCAG8 Antibody

Catalog No: #43294

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	SDCCAG8 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total SDCCAG8 protein.
Immunogen Type	peptide
Immunogen Description	Synthetic peptide of human SDCCAG8
Target Name	SDCCAG8
Other Names	BBS16; CCCAP; SLSN7; NPHP10; hCCCAP; HSPC085; NY-CO-8; CCCAP SLSN7
Accession No.	Swiss-Prot#: Q86SQ7 Gene ID: 10806
Uniprot	Q86SQ7
GeneID	10806;
Concentration	0.9mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:20-1:100

Images



Immunohistochemical analysis of paraffin-embedded Human esophagus cancer tissue using #43294 at dilution 1/20,

Background

This gene encodes a centrosome associated protein. This protein may be involved in organizing the centrosome during interphase and mitosis. Mutations in this gene are associated with retinal-renal ciliopathy.

Note: This product is for in vitro research use only