

## CDHR4 Antibody

Catalog No: #43364

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	CDHR4 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total CDHR4 protein.
Immunogen Description	Fusion protein of human CDHR4
Target Name	CDHR4
Other Names	CDH29; PRO34300
Accession No.	Swiss-Prot#: A6H8M9Gene ID: 389118
Uniprot	A6H8M9
GeneID	389118;
Concentration	1.5mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol.
Storage	Store at -20°C

## Application Details

Immunohistochemistry: 1:20-1:100

## Images



Immunohistochemical analysis of paraffin-embedded Human brain tissue using #43364 at dilution 1/25,

## Background

Cadherins are calcium-dependent cell adhesion proteins. They preferentially interact with themselves in a homophilic manner in connecting cells; cadherins may thus contribute to the sorting of heterogeneous cell types (By similarity).

---

Note: This product is for in vitro research use only