

Chicken Immunoglobulin G (IgG) ELISA Kit

Catalog No: #EK18148



Package Size: #EK18148-1 48T #EK18148-2 96T

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Chicken Immunoglobulin G (IgG) ELISA Kit
Brief Description	Chicken Immunoglobulin G (IgG) ELISA Kit
Applications	ELISA
Species Reactivity	Chicken

Application Details

Detect Range:0.625 ug/mL - 20 ug/mL

Sample Type: Serum, plasma, Urine,tissue homogenates and other biological fluids.

Sample Volume:50ul

Assay Time: 3-5h

Detection wavelength:450nm

Note: This product is for in vitro research use only