

Goat Anti-Mouse IgG Secondary Antibody AP Conjugated

Catalog No: #L3033



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Goat Anti-Mouse IgG Secondary Antibody AP Conjugated
Clonality	Polyclonal
Purification	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Species Reactivity	Ms
Specificity	IgG (H+L)
Immunogen Description	Mouse IgG
Conjugates	Alkaline Phosphatase
Formulation	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide
Storage	Store at -20 °C/1 year

Application Details

Antibody Format:Whole IgG
ELISA: 1:5,000-1:50,000 Western Blot: 1:5,000-1:50,000

Note: This product is for in vitro research use only