

TGM6 Antibody

Catalog No: #43425

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TGM6 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total TGM6 protein.
Immunogen Description	Synthetic peptide of human TGM6
Target Name	TGM6
Other Names	TG6; TGY; SCA35; TGM3L; dJ734P14.3
Accession No.	Swiss-Prot#: O95932Gene ID: 343641
Uniprot	O95932
GeneID	343641;
Concentration	2.1mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:20-1:100

Images



Immunohistochemical analysis of paraffin-embedded Human lung cancer tissue using #43425 at dilution 1/35,

Background

The protein encoded by this gene belongs to the transglutaminase superfamily. It catalyzes the cross-linking of proteins and the conjugation of polyamines to proteins. Mutations in this gene are associated with spinocerebellar ataxia type 35 (SCA35). Alternatively spliced transcript variants encoding different isoforms have been found for this gene.?

Note: This product is for in vitro research use only