

GRK2 (Phospho-Ser685) Polyclonal Antibody

Catalog No: #12397

Package Size: #12397-1 50ul #12397-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	GRK2 (Phospho-Ser685) Polyclonal Antibody
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous GRK2 (Phospho Ser685) protein.
Immunogen Type	Peptide
Immunogen Description	Phospho peptide at Ser685 of GRK2
Modification	Phospho
Other Names	ADRBK1; BARK; BARK1; GRK2; Beta-adrenergic receptor kinase 1; Beta-ARK-1
Accession No.	Swiss-Prot#: P25098NCBI Gene ID: 156
Uniprot	P25098
GeneID	156;
Calculated MW	80kd
Concentration	1mg/ml
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Store at -20°C/1 year

Application Details

Western blotting: 1:500-2000

ELISA: 1: 5000

Images



Western Blot analysis of JK cells using Phospho-GRK 2 (S685) Polyclonal Antibody

Note: This product is for in vitro research use only