

DEDD2 Antibody

Catalog No: #43513

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	DEDD2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total DEDD2 protein.
Immunogen Type	protein
Immunogen Description	Fusion protein of human DEDD2
Target Name	DEDD2
Other Names	FLAME-3
Accession No.	Swiss-Prot#: Q8WXF8NCBI Gene ID: 162989
Uniprot	Q8WXF8
GeneID	162989;
Concentration	0.6mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C


Application Details

Immunohistochemistry: 1: 20-100

Images



The image on the left is immunohistochemistry of paraffin-embedded Human colorectal cancer tissue using DEDD2 Antibody at dilution 1/30, on the right is treated with fusion protein. (Original magnification: x200)



The image on the left is immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using DEDD2 Antibody at dilution 1/30, on the right is treated with fusion protein. (Original magnification: x200)

Background

This gene encodes a nuclear-localized protein containing a death effector domain (DED). The encoded protein may regulate the trafficking of caspases and other proteins into the nucleus during death receptor-induced apoptosis. Alternative splicing results in multiple transcript variants.

Note: This product is for in vitro research use only