

ZMAT2 Antibody

Catalog No: #43571



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description


Product Name	ZMAT2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC WB
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total ZMAT2 protein.
Immunogen Type	protein
Immunogen Description	Full length fusion protein
Target Name	ZMAT2
Other Names	Snu23; Ptg-12; hSNU23
Accession No.	Swiss-Prot#: Q96NC0NCBI Gene ID: 153527
Uniprot	Q96NC0
GeneID	153527;
Calculated MW	24kd
Concentration	0.5mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN3, 40% Glycerol.
Storage	Store at -20°C

Application Details


Western blotting: 1:200-1000

Immunohistochemistry: 1: 30-150


Images



Gel: 12%SDS-PAGE
Lysate: 40 µg, Lane 1-6: Jurkat,231,Hela,k562 and A172 cell,Mouse brain tissue lysates,
Primary antibody:ZMAT2 antibody at dilution 1/200,
Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution,
Exposure time: 3 seconds



The image on the left is immunohistochemistry of paraffin-embedded Human tonsil tissue using ZMAT2 Antibody at dilution 1/40, on the right is treated with fusion protein. (Original magnification: x200)



The image on the left is immunohistochemistry of paraffin-embedded Human breast cancer tissue using ZMAT2 Antibody at dilution 1/40, on the right is treated with fusion protein. (Original magnification: x200)

Background

The function of ZMAT2 remains unknown.

Note: This product is for in vitro research use only