

UPK3A Antibody

Catalog No: #43619

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	UPK3A Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total UPK3A protein.
Immunogen Type	protein
Immunogen Description	Fusion protein of human UPK3A
Target Name	UPK3A
Other Names	UP3A; UPK3; UPIII; UPIIIA
Accession No.	Swiss-Prot#: O75631NCBI Gene ID: 7380
Uniprot	O75631
GeneID	7380;
Concentration	1.3mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1: 100-200

Images



The image on the left is immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using UPK3A Antibody at dilution 1/100, on the right is treated with fusion protein. (Original magnification: x200)

Background

This gene encodes a member of the uroplakin family, a group of transmembrane proteins that form complexes on the apical surface of the bladder epithelium. Mutations in this gene may be associated with renal dysplasia. Alternatively spliced transcript variants have been described.

Note: This product is for in vitro research use only