

UTP3 Antibody

Catalog No: #43621

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	UTP3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total UTP3 protein.
Immunogen Type	protein
Immunogen Description	Fusion protein of human UTP3
Target Name	UTP3
Other Names	CRL1; CRLZ1; SAS10
Accession No.	Swiss-Prot#: Q9NQZ2NCBI Gene ID: 57050
Uniprot	Q9NQZ2
GeneID	57050;
Concentration	1.2mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C


Application Details

Immunohistochemistry: 1: 100-200

Images



The image on the left is immunohistochemistry of paraffin-embedded Human colorectal cancer tissue using UTP3 Antibody at dilution 1/100, on the right is treated with fusion protein. (Original magnification: x200)



The image on the left is immunohistochemistry of paraffin-embedded Human tonsil tissue using UTP3 Antibody at dilution 1/100, on the right is treated with fusion protein. (Original magnification: x200)

Background

Essential for gene silencing: has a role in the structure of silenced chromatin. Plays a role in the developing brain (By similarity).

Note: This product is for in vitro research use only