

CIZ1 Conjugated Antibody

Catalog No: #C37856



Package Size: #C37856-AF350 100ul #C37856-AF405 100ul #C37856-AF488 100ul

#C37856-AF555 100ul #C37856-AF594 100ul #C37856-AF647 100ul

#C37856-AF680 100ul #C37856-AF750 100ul #C37856-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CIZ1 Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total CIZ1 protein.
Immunogen Description	Synthetic peptide corresponding to a region derived from internal residues of human CDKN1A interacting zinc finger protein 1
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	NP94; LSFR1; ZNF356
Accession No.	Swiss-Prot#:Q9ULV3NCBI Gene ID:25792NCBI mRNA#:NCBI Protein#:NP_115520/Q9H0N0
Uniprot	Q9ULV3
GeneID	25792;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	100
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Background

The protein encoded by this gene is a zinc finger DNA binding protein that interacts with CIP1, part of a complex with cyclin E. The encoded protein may regulate the cellular localization of CIP1. Several transcript variants encoding different isoforms have been found for this gene.?

Note: This product is for in vitro research use only