

Gzmd Antibody

Catalog No: #43693

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Gzmd Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	WB
Species Reactivity	Ms
Specificity	The antibody detects endogenous levels of total Gzmd protein.
Immunogen Type	peptide
Immunogen Description	Synthetic peptide of mouse Gzmd
Target Name	Gzmd
Other Names	CCP5; Ctla5; Ctla-5
Accession No.	Swiss-Prot#: P11033 NCBI Gene ID: 14941
Uniprot	P11033
GeneID	14941;
Calculated MW	28kd
Concentration	0.4mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:200-1000

Images



Gel: 12%SDS-PAGE

Lysate: 40 µg, Lane: Mouse heart tissue,

Primary antibody: Gzmd antibody at dilution 1/350 dilution,

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution,

Exposure time: 10 seconds

Background

This enzyme is probably necessary for target cell lysis in cell-mediated immune responses.

Note: This product is for in vitro research use only