

APOBEC3A Conjugated Antibody

Catalog No: #C44289



Package Size: #C44289-AF350 100ul #C44289-AF405 100ul #C44289-AF488 100ul

#C44289-AF555 100ul #C44289-AF594 100ul #C44289-AF647 100ul

#C44289-AF680 100ul #C44289-AF750 100ul #C44289-Biotin 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	APOBEC3A Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Species Reactivity	Hu Ms
Specificity	APOBEC3A Antibody detects endogenous levels of total APOBEC3A
Immunogen Description	A synthesized peptide derived from human APOBEC3A
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	APOBEC3A, ARP3, BK150C2.1, Phorbolin i, Phorbolin-1, PHRBN, Phorbolin 1
Accession No.	Swiss-Prot#:P31941NCBI Gene ID:100913187/200315NCBI mRNA#:NCBI Protein#:
Uniprot	P31941
GeneID	100913187;200315;
Excitation Emission	AF350: 346nm/442nm
	AF405: 401nm/421nm
	AF488: 493nm/519nm
	AF555: 555nm/565nm
	AF594: 591nm/614nm
	AF647: 651nm/667nm
	AF680: 679nm/702nm
	AF750: 749nm/775nm
Calculated MW	26
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only