

CACNG8 Conjugated Antibody

Catalog No: #C44701



Package Size: #C44701-AF350 100ul #C44701-AF405 100ul #C44701-AF488 100ul

#C44701-AF555 100ul #C44701-AF594 100ul #C44701-AF647 100ul

#C44701-AF680 100ul #C44701-AF750 100ul #C44701-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CACNG8 Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Species Reactivity	Hu Ms Rt
Specificity	CACNG8 antibody detects endogenous levels of total CACNG8
Immunogen Description	A synthesized peptide derived from human CACNG8
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	CACNG antibody CACNG6 antibody Cacng8 antibody Calcium channel voltage dependent gamma subunit 8 antibody CCG8_HUMAN antibody Neuronal voltage gated calcium channel gamma 8 subunit antibody Neuronal voltage-gated calcium channel gamma-8 subunit antibody TARP gamma-8 antibody Transmembrane AMPAR regulatory protein gamma-8 antibody Voltage dependent calcium channel gamma 8 subunit antibody Voltage-dependent calcium channel gamma-8 subunit antibody
Accession No.	Swiss-Prot#:Q8WXS5NCBI Gene ID:59283NCBI mRNA#:NCBI Protein#:
Uniprot	Q8WXS5
GeneID	59283;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	43
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only