

CITED2 Conjugated Antibody

Catalog No: #C44771

Package Size: #C44771-AF350 100ul #C44771-AF405 100ul #C44771-AF488 100ul

#C44771-AF555 100ul #C44771-AF594 100ul #C44771-AF647 100ul

#C44771-AF680 100ul #C44771-AF750 100ul #C44771-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CITED2 Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Species Reactivity	Hu Ms Rt
Specificity	CITED2antibody detects endogenous levels of total?CITED2
Immunogen Description	A synthesized peptide derived from human CITED2
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	ASD8 antibody Cbp/p300 interacting transactivator 2 antibody Cbp/p300 interacting transactivator with Glu/Asp rich carboxy terminal domain 2 antibody Cbp/p300-interacting transactivator 2 antibody CITE2_HUMAN antibody CITED 2 antibody CITED2 antibody melanocyte-specific gene 1-related gene 1 antibody MRG 1 antibody MRG-1 antibody MRG1 antibody MSG related protein 1 antibody MSG-related protein 1 antibody MSG1-related gene 1 antibody P35srj antibody VSD2 antibody
Accession No.	Swiss-Prot#:Q99967NCBI Gene ID:10370NCBI mRNA#:NCBI Protein#:
Uniprot	Q99967
GeneID	10370;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	28

Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only