

OR4C12 Conjugated Antibody

Catalog No: #C45442



Package Size: #C45442-AF350 100ul #C45442-AF405 100ul #C45442-AF488 100ul
 #C45442-AF555 100ul #C45442-AF594 100ul #C45442-AF647 100ul
 #C45442-AF680 100ul #C45442-AF750 100ul #C45442-Biotin 100ul

Orders: order@signalwayantibody.com
 Support: tech@signalwayantibody.com

Description

Product Name	OR4C12 Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Species Reactivity	Hu
Specificity	OR4C12 Antibody detects endogenous levels of total OR4C12
Immunogen Description	A synthesized peptide derived from human OR4C12
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	OR4C12, Olfactory receptor OR11-259, OR11-259, Olfactory receptor 4C12
Accession No.	Swiss-Prot#:Q96R67NCBI Gene ID:283093NCBI mRNA#:NCBI Protein#:
Calculated MW	34
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.