

CT47A Conjugated Antibody

Catalog No: #C45530



Package Size: #C45530-AF350 100ul #C45530-AF405 100ul #C45530-AF488 100ul

#C45530-AF555 100ul #C45530-AF594 100ul #C45530-AF647 100ul

#C45530-AF680 100ul #C45530-AF750 100ul #C45530-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CT47A Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Species Reactivity	Hu
Specificity	CT47A Antibody detects endogenous levels of total CT47A
Immunogen Description	A synthesized peptide derived from human CT47A
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Cancer/testis antigen 47A Cancer/testis antigen 47 CT47
Accession No.	Swiss-Prot#:Q5JQC4NCBI Gene ID:100507170/255313/653282/728036/728042/728049/728062/728072/728075/728082/728090/728096NCBI mRNA#:NCBI Protein#:
Uniprot	Q5JQC4
GeneID	100507170;255313;653282;728036;728042;728049;728062;728072;728075;728082;728090;728096;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	30
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only