

## CNOT4 Conjugated Antibody

Catalog No: #C46524



Package Size: #C46524-AF350 100ul #C46524-AF405 100ul #C46524-AF488 100ul

#C46524-AF555 100ul #C46524-AF594 100ul #C46524-AF647 100ul

#C46524-AF680 100ul #C46524-AF750 100ul #C46524-Biotin 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |                                                                                                                                                                              |
|-----------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | CNOT4 Conjugated Antibody                                                                                                                                                    |
| Host Species          | Rabbit                                                                                                                                                                       |
| Clonality             | Polyclonal                                                                                                                                                                   |
| Species Reactivity    | Hu Ms                                                                                                                                                                        |
| Specificity           | The antibody detects endogenous levels of total CNOT4 protein.                                                                                                               |
| Immunogen Description | Synthetic peptide corresponding to residues near the N terminal of human CNOT4                                                                                               |
| Conjugates            | Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750                                                                                                                       |
| Other Names           | NOT4; NOT4H; CLONE243                                                                                                                                                        |
| Accession No.         | Swiss-Prot#:O95628NCBI Gene ID:4850NCBI mRNA#:NCBI Protein#:NP_001177777                                                                                                     |
| Uniprot               | O95628                                                                                                                                                                       |
| GeneID                | 4850;                                                                                                                                                                        |
| Excitation Emission   | AF350: 346nm/442nm<br>AF405: 401nm/421nm<br>AF488: 493nm/519nm<br>AF555: 555nm/565nm<br>AF594: 591nm/614nm<br>AF647: 651nm/667nm<br>AF680: 679nm/702nm<br>AF750: 749nm/775nm |
| Calculated MW         | 64                                                                                                                                                                           |
| Formulation           | 0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide                                                                                  |
| Storage               | Store at 4°C in dark for 6 months                                                                                                                                            |

## Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

## Background

---

The protein encoded by this gene is a subunit of the CCR4-NOT complex, a global transcriptional regulator. The encoded protein interacts with CNOT1 and has E3 ubiquitin ligase activity. Several transcript variants encoding different isoforms have been found for this gene.

---

Note: This product is for in vitro research use only