

## URGCP Antibody

Catalog No: #43989

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	URGCP Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total URGCP protein.
Immunogen Type	peptide
Immunogen Description	Synthetic peptide of human URGCP
Target Name	URGCP
Other Names	URG4
Accession No.	Swiss-Prot#: Q8TCY9NCBI Gene ID: 55665
Uniprot	Q8TCY9
GeneID	55665;
Concentration	0.5mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol.
Storage	Store at -20°C

## Application Details

Immunohistochemistry: 1: 20-100

## Images



The image on the left is immunohistochemistry of paraffin-embedded Human liver cancer tissue using URGCP Antibody at dilution 1/25, on the right is treated with synthetic peptide. (Original magnification: x200)

## Background

URG4 is upregulated in the presence of hepatitis B virus (HBV)-encoded X antigen (HBxAg) and may contribute to the development of hepatocellular carcinoma by promoting hepatocellular growth and survival (Tufan et al., 2002 [PubMed 12082552]).

---

Note: This product is for in vitro research use only