

DYNAP Antibody

Catalog No: #44001

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	DYNAP Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total DYNAP protein.
Immunogen Type	peptide
Immunogen Description	Synthetic peptide of human DYNAP
Target Name	DYNAP
Other Names	C18orf26
Accession No.	Swiss-Prot#: Q8N1N2NCBI Gene ID: 284254
Uniprot	Q8N1N2
GeneID	284254;
Concentration	1.1mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1: 30-150

Images



The image on the left is immunohistochemistry of paraffin-embedded Human brain tissue using DYNAP Antibody at dilution 1/50, on the right is treated with synthetic peptide. (Original magnification: x200)

Background

Plays a role in the regulation of cell proliferation. Promotes activation of the AKT1 signaling pathway. Promotes phosphorylation of AKT1 at 'Ser-473'.

Note: This product is for in vitro research use only