

S-adenosylmethionine decarboxylase proenzyme 1

Catalog No: #AP79720



Package Size: #AP79720-1 50ug #AP79720-2 100ug #AP79720-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	S-adenosylmethionine decarboxylase proenzyme 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	AdoMetDC 1
Accession No.	Uniprot ID:P0DMN7
Uniprot	P0DMN7
GeneID	11702;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Amd1

Note: This product is for in vitro research use only