

Guanine nucleotide-binding protein G(s) subunit alpha isoforms short

Catalog No: #AP79721

Package Size: #AP79721-1 50ug #AP79721-2 100ug #AP79721-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Guanine nucleotide-binding protein G(s) subunit alpha isoforms short
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Alternative gene product encoded by XL-exon,GNAS1
Accession No.	Uniprot ID:P84996
Uniprot	P84996
GeneID	2778;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:GNAS

Note: This product is for in vitro research use only