

Proton-associated sugar transporter A

Catalog No: #AP79735

Package Size: #AP79735-1 50ug #AP79735-2 100ug #AP79735-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Proton-associated sugar transporter A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	PAST-A,Solute carrier family 45 member 1,Dnb5
Accession No.	Uniprot ID:Q8K4S3
Uniprot	Q8K4S3
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Slc45a1

Note: This product is for in vitro research use only