

G-protein coupled bile acid receptor 1

Catalog No: #AP79765

Package Size: #AP79765-1 50ug #AP79765-2 100ug #AP79765-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	G-protein coupled bile acid receptor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Membrane-type receptor for bile acids,M-BAR,Tgr5
Accession No.	Uniprot ID:Q80T02
Uniprot	Q80T02
GeneID	338443;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Gpbar1

Note: This product is for in vitro research use only