

# Melanoma-derived growth regulatory protein

Catalog No: #AP79767

Package Size: #AP79767-1 50ug #AP79767-2 100ug #AP79767-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Melanoma-derived growth regulatory protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Cartilage-derived retinoic acid-sensitive protein,CD-RAP,Melanoma inhibitory activity protein,CDRAP
Accession No.	Uniprot ID:Q28038
Uniprot	Q28038
GeneID	280857;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:MIA

Note: This product is for in vitro research use only