

Alpha-tocopherol transfer protein-like

Catalog No: #AP79774

Package Size: #AP79774-1 50ug #AP79774-2 100ug #AP79774-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Alpha-tocopherol transfer protein-like
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Accession No.	Uniprot ID:Q9D3D0
Uniprot	Q9D3D0
GeneID	76080;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Ttpal

Note: This product is for in vitro research use only