

Pro-thyrotropin-releasing hormone

Catalog No: #AP79776

Package Size: #AP79776-1 50ug #AP79776-2 100ug #AP79776-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Pro-thyrotropin-releasing hormone
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Pro-TRH, Prothyroliberin
Accession No.	Uniprot ID:Q62361
Uniprot	Q62361
GeneID	22044;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Trh

Note: This product is for in vitro research use only