



Embryonic stem cell-specific  
5-hydroxymethylcytosine-binding protein

Catalog No: #AP79799

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP79799-1 50ug #AP79799-2 100ug #AP79799-3 1mg

Description

Product Name	Embryonic stem cell-specific 5-hydroxymethylcytosine-binding protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	ES cell-specific 5hmC-binding protein,Putative peptidase SRAPD1,SRAP domain-containing protein 1,C3orf37,DC12,SRAPD1
Accession No.	Uniprot ID:Q96FZ2
Uniprot	Q96FZ2
GeneID	56941;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:HMCES

Note: This product is for in vitro research use only