

Fibroblast growth factor 7

Catalog No: #AP79816



Package Size: #AP79816-1 50ug #AP79816-2 100ug #AP79816-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Fibroblast growth factor 7
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	FGF-7,Heparin-binding growth factor 7,HBGF-7,Keratinocyte growth factor,Fgf-7,Kgf
Accession No.	Uniprot ID:Q02195
Uniprot	Q02195
GeneID	29348;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Fgf7

Note: This product is for in vitro research use only