

Fibroblast growth factor 7

Catalog No: #AP79816



Package Size: #AP79816-1 50ug #AP79816-2 100ug #AP79816-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Fibroblast growth factor 7 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Rt |
| Other Names | FGF-7,Heparin-binding growth factor 7,HBGF-7,Keratinocyte growth factor,Fgf-7,Kgf |
| Accession No. | Uniprot ID:Q02195 |
| Uniprot | Q02195 |
| GeneID | 29348; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Fgf7

Note: This product is for in vitro research use only