

Ezrin

Catalog No: #AP79838



Package Size: #AP79838-1 50ug #AP79838-2 100ug #AP79838-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Ezrin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Cytovillin,Villin-2,p81,VIL2
Accession No.	Uniprot ID:P15311
Uniprot	P15311
GeneID	7430;
SDS-PAGE MW	13.75kDa
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Expression Region:291-415aaStorage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:EZR

Note: This product is for in vitro research use only