

Carboxyl-terminal PDZ ligand of neuronal nitric oxide synthase protein

Catalog No: #AP79849

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP79849-1 50ug #AP79849-2 100ug #AP79849-3 1mg

Description

Product Name	Carboxyl-terminal PDZ ligand of neuronal nitric oxide synthase protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	C-terminal PDZ ligand of neuronal nitric oxide synthase protein,Nitric oxide synthase 1 adaptor protein,Capon
Accession No.	Uniprot ID:O54960
Uniprot	O54960
GeneID	192363;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Nos1ap

Note: This product is for in vitro research use only