

Phosphatidylinositol 5-phosphate 4-kinase type-2 gamma

Catalog No: #AP79855

Package Size: #AP79855-1 50ug #AP79855-2 100ug #AP79855-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Phosphatidylinositol 5-phosphate 4-kinase type-2 gamma
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Phosphatidylinositol 5-phosphate 4-kinase type II gamma,PI(5)P 4-kinase type II gamma,Pip5k2c
Accession No.	Uniprot ID:O88370
Uniprot	O88370
GeneID	140607;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Pip4k2c

Note: This product is for in vitro research use only