

Phosphatidylinositol N-acetylglucosaminyltransferase subunit Y

Catalog No: #AP79856

Package Size: #AP79856-1 50ug #AP79856-2 100ug #AP79856-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Phosphatidylinositol N-acetylglucosaminyltransferase subunit Y
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Phosphatidylinositol-glycan biosynthesis class Y protein,PIG-Y,Pigyl
Accession No.	Uniprot ID:P0C1P0
Uniprot	P0C1P0
GeneID	66268;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Pigy

Note: This product is for in vitro research use only