

# Cross-linked C-terminal telopeptide of Collagen alpha-1(I) chain

Catalog No: #AP79878

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP79878-1 50ug #AP79878-2 100ug #AP79878-3 1mg

## Description

Product Name	Cross-linked C-terminal telopeptide of Collagen alpha-1(I) chain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:ICTP

Note: This product is for in vitro research use only