



Phospholysine phosphohistidine inorganic pyrophosphate phosphatase

Catalog No: #AP79900

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP79900-1 50ug #AP79900-2 100ug #AP79900-3 1mg

Description

Product Name	Phospholysine phosphohistidine inorganic pyrophosphate phosphatase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Accession No.	Uniprot ID:Q9D7I5
Uniprot	Q9D7I5
GeneID	76429;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Lhpp

Note: This product is for in vitro research use only