

Protein Lines homolog 1

Catalog No: #AP79901

Package Size: #AP79901-1 50ug #AP79901-2 100ug #AP79901-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Protein Lines homolog 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Protein Lines homolog 2,Wnt-signaling molecule Lines homolog 1,Lins,Lins2,Wins2
Accession No.	Uniprot ID:Q3U1D0
Uniprot	Q3U1D0
GeneID	72635;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Lins1

Note: This product is for in vitro research use only