



Phospholipid hydroperoxide glutathione peroxidase, mitochondrial

Catalog No: #AP79908

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Package Size: #AP79908-1 50ug #AP79908-2 100ug #AP79908-3 1mg

Description

Product Name	Phospholipid hydroperoxide glutathione peroxidase, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Glutathione peroxidase 4,GPx-4
Accession No.	Uniprot ID:Q91XR8
Uniprot	Q91XR8
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Gpx4

Note: This product is for in vitro research use only