

14-3-3 protein sigma

Catalog No: #AP79996



Package Size: #AP79996-1 50ug #AP79996-2 100ug #AP79996-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	14-3-3 protein sigma
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Epithelial cell marker protein 1,Stratifin,HME1
Accession No.	Uniprot ID:P31947
Uniprot	P31947
GeneID	2810;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SFN

Note: This product is for in vitro research use only