

6-phosphogluconate dehydrogenase, decarboxylating

Catalog No: #AP79998



Package Size: #AP79998-1 50ug #AP79998-2 100ug #AP79998-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	6-phosphogluconate dehydrogenase, decarboxylating
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	PGDH
Accession No.	Uniprot ID:P52209
Uniprot	P52209
GeneID	5226;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: PGD

Note: This product is for in vitro research use only