

14-3-3 protein theta

Catalog No: #AP80005

Package Size: #AP80005-1 50ug #AP80005-2 100ug #AP80005-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	14-3-3 protein theta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	14-3-3 protein tau
Accession No.	Uniprot ID:P68255
Uniprot	P68255
GeneID	25577;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Ywhaq

Note: This product is for in vitro research use only