

3-hydroxyisobutyrate dehydrogenase, mitochondrial

Catalog No: #AP80015

Package Size: #AP80015-1 50ug #AP80015-2 100ug #AP80015-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	3-hydroxyisobutyrate dehydrogenase, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	HIBADH
Accession No.	Uniprot ID:Q2HJD7
Uniprot	Q2HJD7
GeneID	512002;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: HIBADH

Note: This product is for in vitro research use only