

Kynurenine/alpha-aminoadipate aminotransferase, mitochondrial



Catalog No: #AP80017

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Package Size: #AP80017-1 50ug #AP80017-2 100ug #AP80017-3 1mg

Description

Product Name	Kynurenine/alpha-aminoadipate aminotransferase, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	KAT/AadAT,2-aminoadipate aminotransferase,2-aminoadipate transaminase,Alpha-aminoadipate aminotransferase,AadAT,Kynurenine aminotransferase II,Kynurenine--oxoglutarate aminotransferase II,Kynurenine--oxoglutarate transaminase 2,Kynurenine--oxoglutarate transaminase II,KAT2
Accession No.	Uniprot ID:Q8N5Z0
Uniprot	Q8N5Z0
GeneID	51166;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:AADAT

Note: This product is for in vitro research use only