

## 4F2 cell-surface antigen heavy chain

Catalog No: #AP80069

Package Size: #AP80069-1 50ug #AP80069-2 100ug #AP80069-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	4F2 cell-surface antigen heavy chain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	4F2hc,Solute carrier family 3 member 2,CD98
Accession No.	Uniprot ID:Q794F9
Uniprot	Q794F9
GeneID	50567;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Slc3a2

Note: This product is for in vitro research use only