

N-terminal propeptide of Collagen alpha-1(III) chain

Catalog No: #AP80087



Package Size: #AP80087-1 50ug #AP80087-2 100ug #AP80087-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	N-terminal propeptide of Collagen alpha-1(III) chain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Accession No.	Uniprot ID:P13941
Uniprot	P13941
GeneID	84032;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PIIINP

Note: This product is for in vitro research use only