



5'-AMP-activated protein kinase catalytic subunit  
alpha-2

Catalog No: #AP80093

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP80093-1 50ug #AP80093-2 100ug #AP80093-3 1mg

Description

Product Name	5'-AMP-activated protein kinase catalytic subunit alpha-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	AMPK subunit alpha-2,Acetyl-CoA carboxylase kinase,ACACA kinase,Hydroxymethylglutaryl-CoA reductase kinase,HMGCR kinase,Ampk,Ampk2
Accession No.	Uniprot ID:Q09137
Uniprot	Q09137
GeneID	78975;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Prkaa2

Note: This product is for in vitro research use only