

Alpha-1-antitrypsin

Catalog No: #AP80096



Package Size: #AP80096-1 50ug #AP80096-2 100ug #AP80096-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Alpha-1-antitrypsin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	Alpha-1 protease inhibitor, Alpha-1-antiproteinase, Serpin A1, PI
Accession No.	Uniprot ID:P50447
Uniprot	P50447
GeneID	397688;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: SERPINA1

Note: This product is for in vitro research use only