

Alpha-1B-glycoprotein

Catalog No: #AP80105



Package Size: #AP80105-1 50ug #AP80105-2 100ug #AP80105-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Alpha-1B-glycoprotein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Alpha-1-B glycoprotein
Accession No.	Uniprot ID:Q19LI2
Uniprot	Q19LI2
GeneID	117586;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: A1bg

Note: This product is for in vitro research use only