

Putative aspartate aminotransferase, cytoplasmic 2

Catalog No: #AP80139



Package Size: #AP80139-1 50ug #AP80139-2 100ug #AP80139-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Putative aspartate aminotransferase, cytoplasmic 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Glutamate oxaloacetate transaminase 1-like protein 1,Transaminase A-like protein 1
Accession No.	Uniprot ID:Q8NHS2
Uniprot	Q8NHS2
GeneID	137362;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:GOT1L1

Note: This product is for in vitro research use only